# & Clinical Studies

Experience

Clinical Book

Articles & Publications

Preview

Marouene Ben Hadj Hassine - Paolo Bucci - Roberta Gasparro - Alessandro Espedito  
Gilberto Sammartino - Enguerran Lyautey - David Mailhes - Pantaleo Giuseppe  
Nuzzolo Paolo - Amato Massimo - Francesco Riccitiello - Bernard Chapotat - Eric  
Schneck - Vincenzo Cerone - Oreste Trosin - Alberto Maltagliati - Andrea Ottonello  
Giulio Raffaghello - Andrea Mascolo - Giovanni Falisi - Massimo Galli - Pedro Vittorini  
Velasquez - Juan Carlos Gallegos Rivera - Roberto Minasi - Alberto De Biase - Carlo  
Di Paolo - Olivier Granjon - Enzo Di Iorio - Marco Berardini - Alain Simonpieri - Stephan  
Fraisier - Gilles Montalbot - Caroline Bolle - Marie-Paule Gustin - Didier Fau - Patrick  
Exbrayat - Georges Boivin - Brigitte Grosogeat - Marco Esposito - Carlo Barausse  
Roberto Pistilli - Vittorio Checchi - Michele Diazi - Maria Rosaria Gatto - Pietro Felice

# SOMMAIRE

## In-Kone® & EVL®

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## In-Kone® & EVL®

### Safe approach in “all-on-four” technique: a case report

Marouene Ben Hadj Hassine, DDS

Paolo Bucci, MD, DDS

Roberta Gasparro, DDS, PhD

Alessandro Espedito Di Lauro, DDS

Gilberto Sammartino, MD, DDS

*Annali di Stomatologia 2014;V(4):142-145*

The “All-on-Four” concept is based on the placement of four implants in the anterior part of fully edentulous jaws to support a provisional, fixed and immediately loaded full-arch prosthesis. Combining tilted and straight implants for supporting fixed prostheses can be considered a viable treatment modality resulting in a more simple and less time consuming procedure, in significantly less morbidity, in decreased financial costs and a more comfortable postsurgical period for the patients. The authors present a case report with mandibular atrophy and left mental foramina on the top of the residual crest.

**Keywords:** edentulous jaw, dental implants, implant placement

Preoperative panoramic radiograph

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## Minimally invasive approaches to optimize block grafting: a case report

Gilberto Sammartino, Giuseppe Pantaleo, Paolo Nuzzolo, Massimo Amato, Francesco Riccitiello.

*J Oral Implantol. 2014 Sep 18*

In narrow alveolar ridge many surgical approaches exist, and the most frequently used is bone grafting. The sub-periosteal tunneling approach is a minimal safety procedure that allows to allocate the graft in a space that is obtained between the soft tissues and the underlying bone, through an access represented by one single incision on the mesial limit of the bone defect.

A 55-year-old caucasian woman with a severe bone atrophy in the maxillary arch was treated with horizontal bone augmentation. A muco-periosteal tunnel with a periosteal elevator through one single incision of moderate dimensions per hemi-arch, insertion of tissue bank ilium bi-cortical bone (FFB) and a consequent complete prosthetic-implant rehabilitation was realized. The tunnel technique provided primary intention closure of the surgical wound, avoiding dehiscences and infections, and reducing the edema and the post-operative discomfort for the patient.

Clinical and histological studies support the use of FFB and it has been suggested that the results of a graft with homologous and with autologous bone can be comparable, even though the healing phase is more critical and longer for the allograft. The combination of homologous grafts with the tunnel technique reduces significantly the surgical trauma and the postoperative discomfort, as well as the risk of exposure of the graft during the healing phase, with the risk of graft failure. In our experience, onlay grafting with frozen bone is a predictable technique for horizontal augmentation.

**Keywords :** tunnel technique, block grafting, bone augmentation, fresh frozen bone, PRF

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
## The Platform Switching Approach to Optimize Split Crest Technique. Case report

Gilberto Sammartino, Vincenzo Cerone, Roberta Gasparro, Francesco Riccitiello, Oreste Trosino

*Case Rep Dent.* 2014; 2014: 850470.

Published online 2014 Aug 6. doi: 10.1155/2014/850470

Implant rehabilitation of edentulous sites with bone atrophy represents a situation in which dental implant placement might be complex or impossible if regeneration and bone augmentation techniques are not used. One of these techniques is ERE, introduced by Tatum and subsequently modified by Scipioni et al.

This procedure is mainly indicated in cases with sufficient bone height but inadequate thickness. Anyway at least 2-3mm of initial crestal width is mandatory to perform this technique.


ERE can be followed by simultaneous implant placement in order to maintain the created space. This space can be filled with autologous/heterologous graft, with biomaterials or leaving the clot stabilized by a membrane, with or without the application of platelets concentrates such as PRP or PRF that seems to accelerate the healing of hard and soft tissues.

ERE surgical approach is a reliable technique but does not prevent the peri-implant bone resorption.

Different studies have evaluated the peri-implant bone resorption after implant positioning with ERE technique. Strietzel et al. report that 6 months after functional loading the marginal bone around the implants was reduced on average by 1 mm. In another study, during the same period of observation, the mean bone loss was 2 mm. Jensen et al. reported an average of resorption of 1.57 mm (mesial side) and 1.42 mm (distal side) during an average 4.2 years.

This case shows the absence of bone resorption and a slight bone apposition above the implant in a split crest technique using platform switching associated with Morse-cone connection.

**Keywords:** edentulous, bone augmentation techniques, dental implants

One-year radiographic follow-up

## Early loading versus immediate loading: case examples

Alberto Maltagliati, Andrea Ottonello, Giulio Raffaghello, Andrea Mascolo

*Implant practice; vol.6 Number 1: 40-43*

“Immediate loading” is a favorite way to maintain appropriate anatomic and topographic proportions between dental arches, to maintain the contact with periodontal tissues, and for the “occlusal/ contact memory”.<sup>1-2</sup> An edentulous area that remains this way for a long period of time may modify and even preclude the possibility of recreating an occlusion that is esthetic, functional, always predictable, and similar to the original.<sup>3</sup> The concept of “primary stability,” torque wrench, osseointegration, and adequate/suitable denture, are the fundamental factors for immediate and long-distance clinical results.<sup>4-5</sup> With immediate loading, the management of peri-implant soft tissues can lead to a long-term predictability with 30-year old case-studies;<sup>6-7</sup> moreover, the waiting time for osseointegration allows for modulation/mediation of some human mistakes not managed in the early loading phases in the 72-hour prosthetic technique.<sup>8</sup> The use of autologous fibrin glue in the same surgery increases the neoangiogenesis and the healing factors’ activity, also obtaining a better predictability regarding soft-tissue management.

**Keywords:** immediate loading, soft tissue


Table 1

	Male	Female	Part. edent	Tot. edent
Early loading	8	10	8	10
Immediate loading	7	6	5	8
Total	15	16	13	18

Table 2

	No disease	Infectioud disease	Heart disease	Diabetes
Early loading	12	2	3	1
Immediate loading	11	0	2	0
Total	23	2	5	1

Table 3

## Use of 3D cartilage scaffolds for the stabilization of implants and bone regeneration with the fit-lock technique

Giovanni Falisi, Massimo Galli, Pedro Vittorini-Velasquez, Juan C. Gallegos-Rivera, Roberto Minasi, Alberto De Biase, Carlo Di Paolo

*Acta. Odontol. Latinamo. Vol.26 N°3;2013: 167-172*

The surgical procedures for implant applications on the lateral upper areas depend on sinus pneumatization and availability of the residual bone.

In these cases, autologous bone grafting remains the gold standard. Nevertheless, because of the morbidity associated to the donor site and the post-surgical complications, several alternative bone substitutes have been introduced, which, however, imply additional costs and show limited osteoinductive properties.

Such limitations can be compensated with new regeneration strategies for biological and mechanical tissue restoration, a subject which has been addressed by tissue engineering in recent years.

The authors present a new therapeutic option for implant application in the upper maxilla with bone availability less than 4 mm by using 3D scaffolds obtained from antigen-free porcine cartilage in the fit-lock technique.

A longitudinal study on 18 consecutive cases was performed, with a 95.2% success rate one year after the implant. The advantages of this new technique are: 1) Functional and anatomical recovery of the maxillary antrum, 2) Immediate application of the implants; 3) Reduction of surgical times; 4) Absence of patient morbidity; 5) Local anesthesia; 6) Use of implants with a diameter > 4 mm.

**Keywords :** bone regeneration, dental implants

## Implant stability evaluation by resonance frequency analysis in the fit lock technique. A clinical study

Giovanni Falisi DDS, PhD  
Massimo Galli, MDS  
Pedro Vittorino Velasquez, MDS  
Juan Carlos Gallegos Rivera, DDS  
Carlo Di Paolo, MDS

*Annali di Stomatologia 2013;IV(2):196-203*

Surgical procedures for the application of implants in the lateral-superior sectors are affected by the availability of the residual bone.

When this condition is lower than 5 mm it is recommended that techniques involving two therapeutic phases, a reconstructive and an applicative one, as reported in the international literature, are adopted. The authors propose here a new method with the potential to apply implants simultaneously with the reconstructive phase.

The aim of this longitudinal retrospective study was to evaluate the stability of implants applied with the fit lock technique in the upper maxillary in us with bone availability lower than 4 mm by measuring resonance frequency at different follow-up periods. The same as urements, carried out on 30 implants, were analysed with specific statistical procedures.

The results indicate that the stability of the implants inserted with the fit lock method increases progressively over time in a statistically significant manner. The stability recorded after one year from the insertion (ISQ T2) is significantly higher than that recorded after six months (ISQ T1), and this is significantly higher than that recorded at the time of implant placement (ISQ T0).

The implant inserted in the maxillary zones with scarce bone availability and applied with this technique showed a similar stability as reported with other techniques.

In light of the results, the authors confirm that the primary stability represents the basic requirement to guarantee a correct healing of the implant and demonstrate that the fit lock technique also allows reaching this condition when bone availability is minimal.

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
## Use of platform-switched implants with internal conical anti-rotational connection in single-tooth implant replacement

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**Materials and methods.** 15 cases of monoedentulism were rehabilitated using Universal In-Kone® dental implants (Tekka, Brignais, France) which had internal conical connection, and platform-switching design. Implants were placed 1.5 mm below the alveolar ridge. Periapical radiographs were taken before and after surgical and prosthetic procedures. One year after the occlusal load the following clinical parameters were evaluated by the authors: marginal bone loss, plaque index, peri-implant probing depth, suppuration on probing, bleeding on probing, mechanical complications.

**Results.** No marginal bone losses were observed. Mechanical or biological complications, or reinfections did not occur in any case.

**Conclusions.** The cases presented show how the use of internal conical anti-rotational connection and platform-switched implants allows implant marginal bone preservation. Biological or mechanical complications, due to bacterial colonization within the gap between fixture and abutment, with this type of connection were also avoided.

**Keywords:** dental implants, internal conical connection, platform switching, implant-supported rehabilitation, monoedentulism

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## Early Periimplant Tissue Healing on 1-Piece Implants With a Concave Transmucosal Design: A Histomorphometric Study in Dogs

Caroline Bolle, Marie-Paule Gustin, Didier Fau, Patrick Exbrayat, Georges Boivin, Brigitte Grosgeat.

*Implant Dent.* 2015 Aug 13.


The purpose of our study was to investigate the early healing phase of marginal bone and soft tissues around unloaded 1-piece implants with a concave transmucosal design, in a dog model.

**METHODS:** Twenty-four 1-piece implants with a concave transmucosal neck were inserted 1 mm subcrestally in the mandibular ridge of 8 beagle dogs. Four animals were sacrificed after 3 and 12 weeks of healing. Histomorphometric analysis was performed to measure the height of the periimplant tissues.


**RESULTS AND DISCUSSION:** The overall height of the periimplant mucosa was, respectively, 2.67 and 2.52 mm, after 3 and 12 weeks. In the connective tissue, a soft tissue O-ring seal was observed in the healing area provided by the transmucosal concavity, after 12 weeks. The location of the first bone-to-implant contact facing the implant shoulder was 0.00 and +0.18 mm, respectively, after 3 and 12 weeks of healing. Some bone apposition occurred on the implant shoulder during the healing.

**CONCLUSION:** Within the limits of the present study, a concave transmucosal design in 1-piece implants was associated with a short vertical value of biological width and promoted a mechanical interlocking of the implant body at the connective tissue and marginal bone levels.

**Keywords:** narrow transmucosal design, biological width, histomorphometry

Marginal bone healing. After 3 weeks.

Marginal bone healing. After 12 weeks.

## Posterior jaws rehabilitated with partial prostheses supported by 4.0 x 4.0 mm or by longer implants: Four-month post-loading data from a randomised controlled trial

Marco Esposito, Carlo Barausse, Roberto Pistilli, Vittorio Checchi, Michele Diazzi, Maria Rosaria Gatto, Pietro Felice  
*Eur J Oral Implantol.* 2015

**PURPOSE:** To evaluate whether 4.0 x 4.0 mm dental implants could be an alternative to implants at least 8.5 mm long, which were placed in posterior jaws in the presence of adequate bone volumes.

**MATERIALS AND METHODS:** One hundred and fifty patients with posterior (premolar and molar areas) mandibles having at least 12.5 mm bone height above the mandibular canal or 11.5 mm bone height below the maxillary sinus, were randomised according to a parallel group design, in order to receive one to three 4.0 mm-long implants or one to three implants which were at least 8.5 mm long, at three centres. All implants had a diameter of 4.0 mm. Implants were loaded after 4 months with definitive screw-retained prostheses. Patients were followed up to 4-month post-loading and outcome measures were prosthesis and implant failures, any complications and peri-implant marginal bone level changes.

**RESULTS:** No patients dropped-out before the 4-month evaluation. Three patients experienced the early failures of one 4.0 mm-long implant each, in comparison to two patients who lost one long implant each (difference in proportion = 0.01; 95% CI -0.06 to 0.09; P = 0.50). Consequently, two prostheses in each group could not be delivered as planned (difference in proportion = 0; 95% CI -0.07 to 0.07; P = 0.69), and one patient from each group is still waiting to have their prostheses delivered. Three short implant patients experienced three complications versus two long implant patients (difference in proportion = 0.01; 95% CI -0.06 to 0.09; P = 0.50). There were no statistically significant differences in prosthesis failures, implant failures and complications. Patients with short implants lost on average 0.38 mm of peri-implant bone at 4 months and patients with long mandibular implants lost 0.42 mm. There were no statistically significant differences in bone level changes up to 4 months between short and long implants (mean difference = 0.04 mm; 95% CI: -0.041 to 0.117; P = 0.274).

**CONCLUSIONS:** Four months after loading, 4.0 x 4.0 mm implants achieved similar results as 8.5 x 4.0 mm-long or longer implants in posterior jaws, however 5 to 10 years post-loading data are necessary before reliable recommendations can be made.

**Keywords:** posterior jaws, short dental implants

## Rehabilitation of the atrophic posterior mandible with short (4-mm) implants : a case report

Roberto Pistilli, Carlo Barausse, Luigi Checchi, Pietro Felice  
*Int J Periodontics Restorative Dent.* 2014 Sep-Oct;34(5):713-8.

This case report describes a successful implant-prosthetic rehabilitation of an atrophic posterior mandible with 4-mm-long implants. The patient refused to undergo any reconstructive surgery, and because the available bone up to the inferior alveolar nerve was only 5 mm or less, the patient received four implants of 4-mm length. Four months after implant placement, a provisional prosthesis was put in place; after another 4 months, this was then in turn replaced with a definitive prosthesis. The use of such short implants allows a fixed prosthetic solution without the need for vertically augmenting the mandibular bone. This procedure considerably reduces intra- and postoperative patient discomfort compared with reconstructive surgery for the placement of longer implants. The follow-up time was 1 year after implant loading.

**Keywords:** short dental implants

## Rehabilitation of an atrophic posterior mandible with 4-mm short implants: a 3 year post-loading case report

Carlo Barausse, Roberto Pistilli, Pietro Felice  
*Implant Dent.* 2015 Aug 13.



**Aim:** We describe a successful implant-prosthetic rehabilitation of an atrophic posterior mandible with 4-mm long implants.

### Materials and Methods:

A 62-year-old systemically healthy male was referred for fixed prosthetic rehabilitation of the right posterior mandible. Clinical and radiographic assessments showed an important posterior mandibular atrophy. Computed Tomography (CT) scans revealed 5 mm mean residual bone height above the mandibular canal. The patient refused to undergo any reconstructive surgeries, but he wanted a fixed rehabilitation, so the authors proposed placing 4-mm short transmucosal implants (TwinKon Universal SA2, Global D, Lyon, France). After local anesthesia a full-thickness crestal incision was made and 4 implants (4 mm in length and 4 mm in diameter) were placed in the mandibular right second premolar and first, second and third molar sites.

**Results:** The post-operative course was uneventful. Four months after implant placement, a temporary prosthesis was put in place, which was then replaced by a definitive prosthesis after another four months. The follow-up time is at three years after implants loading.



**Conclusions:** With the limits of this case report the use of such short implants can allow a fixed prosthetic solution without the need for vertically augmenting of the mandibular bone. This procedure considerably reduces intra and postoperative patient discomfort compared to reconstructive surgery for placing longer standard implants.

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